

Contents

- [Settings](#)
- [Identify all files with runs](#)
- [Process each run file](#)
- [Data info](#)
- [Open Data File that contains analysis report from CHN instrument](#)
- [Find index for each type of run](#)
- [Instrument signals for each element](#)
- [Carbon: Carbon signal data.Area2](#)
- [Run CARBON calibration curves](#)
- [Calibration run QA](#)
- [Plot Carbon Calibration curve](#)
- [Linear regression of sulfCarb as a f\(x\) of C signal -- post QA](#)
- [Data report of carbon and nitrogen per filter \(NOT normalized by volume yet\)](#)
- [QA/QC with Buffalo River NIST reference material](#)
- [Export figures to file](#)

```
% Process CHN runs from Elementar Cube analyzer
```

```
clearvars
close all
tic

% Runs Info

cruise = 'P06_2017'; %enter here cruise or campaign code
```

Settings

```
calQA = 'yes'; % assess quality of cal runs and remove suspect ones

exportFigs = 'no';

LineWidths = 1.6;

newColors = get(groot,'DefaultAxesColorOrder');

scaledMADs = 3.5; % number os scaled Median Absolute deviations for...
                % calibration run to be classified as an outlier

if ispc
    homeStr = 'C:/Users/jchavesc';
else
    homeStr = '~';
end
```

Identify all files with runs

```
path2Data = [homeStr '/Documents/CVO/Biogeochem/POC/P06_Runs/'];

cd(path2Data)

if ispc
    delete('._2*') % removes hidden files created by Excel
else
    system('rm ._2*') % removes hidden files created by Excel
end

runFiles = dir('*.xlsx');
```

```
rm: cannot remove '._2*': No such file or directory
```

```
ans =
```

```
1
```

Process each run file

```
for k = 1:length(runFiles)
```

Data info

```
rundID = num2str(k); % run ID (temp) *****
runFileStr = strsplit(lower(runFiles(k).name), '_fac'); % to tag plots
```

Open Data File that contains analysis report from CHN instrument

your data is located

```
sheetStr = 'samples'; % sheet tab name

% open data and put into numeric, text, and 'raw' outputs

[num,txt,row] = xlsread(runFiles(k).name,sheetStr);

% use data from 'raw' import from XLS file
% and remove header row

data = raw(2:end,:);

% extract 1st line of headers
headers = txt(1,1:end);
headers = string(headers); % convert to a vector of strings

% Column for run 'Name'

namesCol = headers == "Name";
```

Find index for each type of run

```
runLength = length(data);

% Preallocate vectors for run class indeces

sulfDex = zeros(size(runLength)); % index of sulfanilamide standards
airDex = zeros(size(runLength)); % index of air blanks
tinCalDex = zeros(size(runLength)); % index tin boats blanks
% (i.e. for running weighed standards)
tinSampDex = zeros(size(runLength)); % index tin circle blanks (i.e. for running filter with samples or filter blanks)
buffRivDex = zeros(size(runLength)); % index Buffalo River Sediment reference runs
sampleDex = zeros(size(runLength)); % index of sample runs
sampleSeq = zeros(size(runLength)); % Sample ID

% Preallocate for tin blanks means from all runs to fill runs
% without blank runs

if k == 1
    allCarbCalBlanks = zeros(k,1);
    allCarbSampBlanks = zeros(k,1);
    allNitroCalBlanks = zeros(k,1);
    allNitroSampBlanks = zeros(k,1);
end

% Column with names/ID entered for each run

runNames = data(:,namesCol);

for i = 1:runLength

    if isnumeric(runNames{i}) % if runNames{i} convert directly to number is a sample ID
```

```

sampleDex(i) = 1;
sampleSeq(i) = runNames{i};

else

    runStr = runNames{i}; % string entered for individual run

    % Cleanout runStr for spaces and make it all lowercase

    runStr = lower(runStr(~isspace(runStr)));

    % clean NaNs

    if isnan(runStr)
        runStr = 'empty';
    end

    % sulfanilamide
    sulfDex(i) = any(regexpi(runStr,'sulfa')); % sulfanilamide standards i.e., any runStr containing 'sulfa...'
    sulfDex = logical(sulfDex);

    % air blanks
    airDex(i) = any(regexpi(runStr,'air')); % air blanks i.e., any runStr containing 'air...'
    airDex = logical(airDex);

    % acidified filter blanks
    % acidFiltDex(i) = any(regexpi(runStr,'acid')); % sulfanilamide standards i.e., any runStr containing 'sulfa...'
    % acidFiltDex = logical(acidFiltDex);

    % Tin blanks (boats and sheets)

    if any(regexpi(runStr,'tin'))

        if any(regexpi(runStr,'35x35')) || any(regexpi(runStr,'30x30'))
            tinSampDex(i) = 1;
        else
            tinCalDex(i) = 1;
        end
    end

    tinSampDex = logical(tinSampDex);
    tinCalDex = logical(tinCalDex);

    % Sample runs numbers stored as string

    if ~isnan(str2double(runStr)) % if run string converts to number

        sampleDex(i) = 1;
        sampleSeq(i) = str2double(runStr);

    end

    % Buffalo River NIST reference

    buffRivDex(i) = any(regexpi(runStr,'buff')); % Buffalo River reference i.e., any runStr containing 'buff...'
    buffRivDex = logical(buffRivDex);

end

sampleDex = logical(sampleDex);
end

```

Instrument signals for each element

Carbon: Carbon signal data.Area2

```

% All carbon areas

area2Col = headers == "Area2";

carbonAreas = data(:,area2Col); % Carbon signals (i.e. Area2)

carbonAreas = cell2num(carbonAreas);

% calibration (tin boats) blanks

```

```

calCarbBlankReps = carbonAreas(tinCalDex); % C blanks for calibration runs

calCarbBlank = mean(calCarbBlankReps); % Mean C blanks for calibration runs ONLY

allCarbCalBlanks(k) = calCarbBlank;

if isnan(calCarbBlank)
    calCarbBlank = nanmean(allCarbCalBlanks(1:k));
end

% Sample tin circle blanks

sampleCarbBlankReps = carbonAreas(tinSampDex); % C tin blanks for sample runs

sampleCarbBlank = nanmean(sampleCarbBlankReps); % Mean C tin blanks for sample runs

allCarbSampBlanks(k) = sampleCarbBlank;

if isnan(sampleCarbBlank)

    sampleCarbBlank = nanmean(allCarbSampBlanks(1:k));

end

```

Run CARBON calibration curves

```

% Carbon -- Amount of C in mg for sulfanilamide runs

weightCol = headers == "WeightVol";

weights = data(:,weightCol);

weights = cell2num(weights);

sulfCarbWeight = (41.8/100).* weights(sulfDex); % mg C per sulfanilamide run

sulfCarbAreas = carbonAreas(sulfDex) - calCarbBlank; % Sulfanilamide areas, blank corrected

if isnan(calCarbBlank) || isnan(sampleCarbBlank)
    stop
end

x = sulfCarbAreas;
y = sulfCarbWeight;

```

Calibration run QA

```

if strcmpi(calQA,'yes')

    areaCarbRatio = x./y; % area to carbon ratio in sulfanilamide stds
    medianCarbRatio = nanmedian(areaCarbRatio);
    stdCarbRatio = std(areaCarbRatio);
    absCarbResid = abs(areaCarbRatio - medianCarbRatio);
    maxCarbResid = max(absCarbResid);
    MADCarbRatio = mad(areaCarbRatio,1);
    modZScores = 0.6745 * ((areaCarbRatio - medianCarbRatio) ./ MADCarbRatio);

    % cal QA plot
    % Tweaks so that figures renders the same in Mac & Linux

    screens = handle(0); %
    mainScreen = screens.MonitorPositions(1,3:4);

    if ismac
        h0i = figure('Position',[mainScreen(1) * 1.06 378 364 636]);
    elseif isunix
        h0i = figure('Position',[mainScreen(1) * 1.06 209 500 900]);
        h0i.Renderer = 'OpenGL';
    end

    subplot(2,1,1)

```

```

L1 = plot(areaCarbRatio,'o:');
hold on

% axes properties

h1 = handle(gca);
h1.XLim = [0 ceil(numel(areaCarbRatio) * 1.06)];
h1.LineWidth = LineWidths;
grid on

h1.TickDir = 'out';
h1.TickLength = h1.TickLength .* 1.5;

L2 = hline(medianCarbRatio);
L2.LineStyle = '-.';
L2.Color = newColors(2,:);

% An outlier is a value that is more than three SCALED median absolute ...
% deviations (MAD) away from the median:
%
% For a random variable vector A made up of N scalar observations,
% the median absolute deviation (MAD) is defined as

%     MAD' = median(abs(A-median(A)));
%
% The scaled MAD is defined as c * MAD
% where c = 1.4826 and is given by:
% -1/(sqrt(2)*erfcinv(3/2)).
%
% [***NOTE***] In Matlab mad(x) gives the MEAN absolute deviation from
% the media & while in R the function mad(x) gives the MEDIAN absolute
% deviation from the median. To get the scaled MAD in Matlab use the
% syntax mad(x,1)

% Graphical implementation of the above with the threshold values

X = [h1.XLim fliplr(h1.XLim)];

c = -1/(sqrt(2)*erfcinv(3/2)); %

M_i = scaledMADs * c * MADCarbRatio;

    y1 = medianCarbRatio - M_i;
    y2 = fliplr(medianCarbRatio + M_i);
Y = [y1 y1 y2 y2];
L3 = patch(X,Y,L2.Color);
L3.FaceAlpha = 0.3;
L3.EdgeColor = 'none';

% Y Lims

f = 4;
h1.YLim = [(medianCarbRatio - (M_i * f)) (medianCarbRatio + (M_i * f))];

%     L1 = plot(areaCarbRatio,'o:');

%hline([meanCarbRatio meanCarbRatio-stdCarbRatio meanCarbRatio+stdCarbRatio],{'b', 'r', 'r'},{'mean+/-s.d.',' ',' '})
xLab1 = xlabel('Run Order');
xLab1.Interpreter = 'latex';
yLab1 = ylabel('Signal:C ratio,  $\frac{S_{\rm C}}{\rm mg\ C}$ ');
yLab1.Interpreter = 'latex';
%title(['C calibration Q/A ' rundID ' run'])

% ID outlier Sulfanilamide runs (Mean )areaCarbRatio

goodCalRuns = abs(modZScores) < scaledMADs; % ~isoutlier(areaCarbRatio,'median',3); % index of cal runs to keep
badCalRuns = abs(modZScores) > scaledMADs; % isoutlier(areaCarbRatio,'median',3); % index of cal runs to keep

xx = 1:numel(areaCarbRatio);

L11 = plot(xx(badCalRuns),areaCarbRatio(badCalRuns),'o:');
if ~isempty(L11)
    L11.Color = newColors(2,:);
end

```

```

x = x(goodCalRuns);
y = y(goodCalRuns);

% Legends

LegL3 = patch(h1.XLim(2) * [0.03 0.03 0.1 0.1],h1.YLim(1) + (range(h1.YLim) * [0.76 0.87 0.87 0.76]),L2.Color);
LegL3.FaceAlpha = 0.3;
LegL3.EdgeColor = 'none';
LegL4 = plot([h1.XLim(2) h1.XLim(2)] .* [0.03 0.1],h1.YLim(1) + (repmat(range(h1.YLim),1,2) .* [0.815 0.815]));
LegL4.LineWidth = LineWidths;
LegL4.Color = L2.Color;
ttL1 = text(h1.XLim(2) * 0.1167,h1.YLim(1) + (range(h1.YLim) * 0.815),'$x_{c-} < \hat{x} < x_{c+}$');
ttL1.Interpreter = 'Latex';
ttL1.FontSize = 18;

tt01 = text(h1.XLim(2) * 0.8,h1.YLim(1) + range(h1.YLim * 0.8),'a');
tt01.Units = 'Inches';
tt01.Position = tt01.Position;
tt01.FontSize = 20;

% Box
bx = patch(h1.XLim(2) * [0.018 0.018 0.55 0.55],h1.YLim(1) + (range(h1.YLim) * [0.73 0.9 0.9 0.73]),'w');
bx.LineWidth = LineWidths;

%Sort Graphic objects


%   h1.SortMethod = 'childorder';
%
%   hP = handle(h1.Children(8));
%
%   uistack(hP,'bottom');

ti01 = title(runFileStr{1});
ti01.FontWeight = 'normal';
ti01.Interpreter = 'none';
ti01.FontSize = 12;


end

```


Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.



Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.



Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.



Plot Carbon Calibration curve

```
% C regression statistics
X=[ones(size(x,1),1) x];
[b,bint,r,rint,stats] = regress(y,X);
```

Linear regression of sulfCarb as a f(x) of C signal -- post QA

```
xx = union(0,x); % linear fit evaluation range
carbModel = flipud(b)'; % polyfit(x,y,1); % Carbon mass vs C signal linear fit
%carbModel(2) = 0; % subtraction of intercept
z = polyval(carbModel,xx);
```

```
subplot(2,1,2);
L33 = plot(x,y,'o');
```

```

L33.Color = newColors(1,:);

grid on
hold on

L4 = plot(sulfCarbAreas(goodCalRuns == 0),sulfCarbWeight(goodCalRuns == 0),'o');

if ~isempty(L4)
    L4.Color = newColors(2,:);
end

if sum(badCalRuns) > 0
    l0 = legend(L4,'Excluded','AutoUpdate','off');
    l0.Location = 'northwest';
end

aa = get(gca,'xlim');
bb = get(gca,'ylim');

L34 = plot(xx,z,'-');
L34.Color = newColors(2,:);

yL2 = ylabel('Carbon, mg');
yL2.Interpreter = 'Latex';

xL2 = xlabel('Carbon Signal,  $S_{\rm C}$ ');
xL2.Interpreter = 'Latex';

% Display R^2 & regression eq on cal curve plot

tt01 = text(0.4 * aa(2),0.22 * bb(2),...
    ['$R^2 = ' num2str(stats(1),4) '$'],...
    'fontSize',16,...
    'Interpreter','latex');
tt01.Units = 'Inches';
tt01.Position(1) = 1;
tt01.HorizontalAlignment = 'left';

newStr = regexp(num2str(b(2)),'e-','split');

m = str2double(newStr{1});
m = round(m,3);

newStr{1} = num2str(m);

regString = ['$ y = ' newStr{1} '\times 10^{' newStr{2}(2) '}x + ' num2str(b(1),2) '$'];

tt02 = text(1.7,0.17,...
    regString,...
    'fontSize',16,...
    'Interpreter','latex');
tt02.Units = 'Inches';
tt02.Position(1) = 1.0;
tt02.HorizontalAlignment = 'left';

h2 = handle(gca);
h2.LineWidth = LineWidths;
h2.TickDir = 'out';
h2.TickLength = h2.TickLength .* 1.5;

tt1 = text(h2.XLim(2) * 0.8,h2.YLim(1) + range(h2.YLim * 0.75),'b');
tt1.Units = 'Inches';
%tt1.Position = [3.3571    2.5904    0];
tt1.FontSize = 20;

% Plot position adjustment

%h2.Units = 'Inches';

h2.Position = h2.Position .* [1.25 1.4 0.95 1];



% Carbon

sampleCarbAreas = carbonAreas(sampleDex) - sampleCarbBlank; % Sample Carbon areas blank corrected (minus tin square foil)



```

```
%sampleCarbWeight = polyval(carbModel,sampleCarbAreas); % calculate sample Carbon in mg
% subtracting the intercept of cal curve
sampleCarbWeight = polyval(carbModel,sampleCarbAreas) - carbModel(2); % calculate sample Carbon in mg
```


Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.





Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.



Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.



Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.



Data report of carbon and nitrogen per filter (NOT normalized by volume yet)

QA/QC with Buffalo River NIST reference material

```
% Carbon percent error based on Buffalo River runs
% NIST certified C % content in Buff R. sed. is 3.351%

%Buffalo River Sediment Ref stds

buffRivMeasCarbon = polyval(carbModel,carbonAreas(buffRivDex) - calCarbBlank) - carbModel(2); % measured C in Buff stds

buffRivExpCarbon = 3.351/100 * weights(buffRivDex); % expected C in Buff stds

% tile figures

if k == 1

    screenSize = get(0,'ScreenSize');

    plotWidth = 350;
    plotHeight = 325;
```

```
        nCols = floor(screenSize(3)/plotWidth);
        nRows = floor(screenSize(4)/plotHeight);

    end

    xPos = mod(k - 1,nCols) * plotWidth;
    if k > nCols
        yPos = (screenSize(4) - (67 + ceil(k/nCols) * plotHeight));
    else
        yPos = (screenSize(4) - (ceil(k/nCols) * plotHeight));
    end
    %test(k,1:3) = [k mod(k - 1,nCols) ceil(k/nCols) ];

    f = figure('Position',[xPos yPos plotWidth plotHeight]);
```



```
plot(buffRivExpCarbon,buffRivMeasCarbon,'o')

    grid on
    hold on

a=max(get(gca,'xlim'));
    set(gca,'xlim',[0 a]);
    set(gca,'ylim',[0 a]);
    line([0 a],[0 a],'color','r')


    xlabel('Buff R. sed. exp. Carbon, mg','fontsize',15)
    ylabel('Buff R. sed. meas. Carbon, mg','fontsize',15)

    error = errperf(buffRivExpCarbon,buffRivMeasCarbon,'rmspe');


    text(a * 0.5, a * 0.2,['RMSE%= ' num2str(error,3)],'fontsize',17)

    til = title(['Buff. River: ' runFileStr{1}],...
        'fontsize',15);
    til.Interpreter = 'none';
    til.FontWeight = 'normal';
```


Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.



Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.



Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.



```

%% Export C concentration data to a file along with sample ID

sampleID = str2double(runNames(sampleDex));

% if sampleID are not being converted to numbers try pass them directly

if sum(isnan(sampleID)) == numel(sampleID)
    sampleID = cell2num(runNames(sampleDex));
elseif any(isnan(sampleID))
    a = runNames(sampleDex);
    sampleID(isnan(sampleID)) = cell2num(a(isnan(sampleID)));
end

if k == 22
    stop
end

runNumber = repmat(k,numel(sampleID),1);

M = [runNumber sampleID sampleCarbWeight];

    if k == 1
        dlmwrite('P06carbonWeights.txt',M,...
            'delimiter',' ','roffset',0,'precision',9)

    else

        dlmwrite('P06carbonWeights.txt',M,'-append',...
            'delimiter',' ','roffset',0,'precision', 9);

    end
end

```

```

%% Export calibration curve and other performance metrics to file

runNumberStr = sprintf('runNumber_%d',k);

runData.(runNumberStr).carbModel = carbModel;
runData.(runNumberStr).buffRMSE = error;
runData.(runNumberStr).buffData = [buffRivExpCarbon buffRivMeasCarbon];
runData.(runNumberStr).calCarbBlankReps = calCarbBlankReps;
runData.(runNumberStr).sampleCarbBlankReps = sampleCarbBlankReps;

if k == length(runFiles)
    save('P06runDataStrct.mat','runData')
end

```

```
end
```

Export figures to file

```
ag = findobj; % all graphical objects

nf = ag(2).Number; %max(ag==fix(ag)); % maximal integer value is number of figures

if strcmpi(exportFigs,'yes')

    %cd('/Users/jchavesc/Documents/CVO/Biogeochem/POC/Runs/Figures')
    %cd('/Users/jchavesc/Documents/CVO/Workshops/POC_mfiles')
    %cd(['~/jchavesc/Documents/CVO/Workshops/POC_mfiles']) % Re direct to where

    cd(['~/Documents/CVO/Workshops/Writing/POCLatexVersion/Figures'])

    for i = 1:nf

        print(i,'-dpdf',['Fig' num2str(i) 'CubeRun' rundID '.pdf'])
        print(i,'-dpng','calCurveFigv2_latest.png','-r600')

    end

end
```

```
disp('Process completed')
toc
```

Process completed
Elapsed time is 253.391700 seconds.

Published with MATLAB® R2019b

Contents

- [Settings](#)
- [Identify all files with runs](#)
- [Process each run file](#)
- [Data info](#)
- [Open Data File that contains analysis report from CHN instrument](#)
- [Find index for each type of run](#)
- [Instrument signals for each element](#)
- [Nitrogen: Nitrogen signal data.Area1](#)
- [Run CARBON calibration curves](#)
- [Calibration run QA](#)
- [Plot Nitrogen Calibration curve](#)
- [Linear regression of sulfNitro as a f\(x\) of C signal -- post QA](#)
- [Carbon and Nitrogen in filter samples](#)
- [Export figures to file](#)

```
% Process CHN runs from Elementar Cube analyzer
```

```
clearvars
close all
tic

% Runs Info

cruise = 'P06_2017'; %enter here cruise or campaign code
```

Settings

```
calQA = 'yes'; % assess quality of cal runs and remove suspect ones

exportFigs = 'no';

LineWidths = 1.6;

newColors = get(groot,'DefaultAxesColorOrder');

scaledMADs = 3.5; % number os scaled Median Absolute deviations for...
                % calibration run to be classified as an outlier

if ispc
    homeStr = 'C:/Users/jchavesc';
else
    homeStr = '-';
end
```

Identify all files with runs

```
path2Data = [homeStr '/Documents/CVO/Biogeochem/POC/P06_Runs/'];

cd(path2Data)

if ispc
    delete('._2*') % removes hidden files created by Excel
else
    system('rm ._2*') % removes hidden files created by Excel
end

runFiles = dir('*.xlsx');
```

```
rm: cannot remove '._2*': No such file or directory

ans =

    1
```

Process each run file

```
for k = 1:length(runFiles)
```

Data info

```
    rundID = num2str(k); % run ID (temp) *****

    runFileStr = strsplit(lower(runFiles(k).name),'_fac'); % to tag plots
```

Open Data File that contains analysis report from CHN instrument

your data is located

```
    sheetStr = 'samples'; % sheet tab name
```

```

% open data and put into numeric, text, and 'raw' outputs

[num,txt,raw] = xlsread(runFiles(k).name,sheetStr);

% use data from 'raw' import from XLS file
% and remove header row

data = raw(2:end,:);

% extract 1st line of headers
headers = txt(1,1:end);
headers = string(headers); % convert to a vector of strings

% Column for run 'Name'

namesCol = headers == "Name";

```

Find index for each type of run

```

runLength = length(data);

% Preallocate vectors for run class indices

sulfDex = zeros(size(runLength)); % index of sulfanilamide standards
airDex = zeros(size(runLength)); % index of air blanks
tinCalDex = zeros(size(runLength)); % index tin boats blanks
% (i.e. for running weighed standards)
tinSampDex = zeros(size(runLength)); % index tin circle blanks (i.e. for running filter with samples or filter blanks)
buffRivDex = zeros(size(runLength)); % index Buffalo River Sediment reference runs
sampleDex = zeros(size(runLength)); % index of sample runs
sampleSeq = zeros(size(runLength)); % Sample ID

% Preallocate for tin blanks means from all runs to fill runs
% without blank runs

if k == 1
    allNitroCalBlanks = zeros(k,1);
    allNitroSampBlanks = zeros(k,1);
end

% Column with names/ID entered for each run

runNames = data(:,namesCol);

for i = 1:runLength

    if isnumeric(runNames{i}) % if runNames{i} convert directly to number is a sample ID

        sampleDex(i) = 1;
        sampleSeq(i) = runNames{i};

    else

        runStr = runNames{i}; % string entered for individual run

        % Cleanout runStr for spaces and make it all lowercase

        runStr = lower(runStr(~isspace(runStr)));

        % clean NaNs

        if isnan(runStr)
            runStr = 'empty';
        end

        % sulfanilamide
        sulfDex(i) = any(regexpi(runStr,'sulfa')); % sulfanilamide standards i.e., any runStr containing 'sulfa...'
        sulfDex = logical(sulfDex);

        % air blanks
        airDex(i) = any(regexpi(runStr,'air')); % air blanks i.e., any runStr containing 'air...'
        airDex = logical(airDex);

        % acidified filter blanks
        % acidFiltDex(i) = any(regexpi(runStr,'acid')); % sulfanilamide standards i.e., any runStr containing 'sulfa...'
        % acidFiltDex = logical(acidFiltDex);

        % Tin blanks (boats and sheets)

        if any(regexpi(runStr,'tin'))

            if any(regexpi(runStr,'35x35')) || any(regexpi(runStr,'30x30'))
                tinSampDex(i) = 1;
            else
                tinCalDex(i) = 1;
            end
        end

        tinSampDex = logical(tinSampDex);
        tinCalDex = logical(tinCalDex);

        % Sample runs numbers stored as string

        if ~isnan(str2double(runStr)) % if run string converts to number

            sampleDex(i) = 1;

```

```

        sampleSeq(i) = str2double(runStr);

    end

    % Buffalo River NIST reference

    buffRivDex(i) = any(regexpi(runStr,'buff')); % Buffalo River reference i.e., any runStr containing 'buff...'
    buffRivDex = logical(buffRivDex);

end

sampleDex = logical(sampleDex);
end

```

Instrument signals for each element

Nitrogen: Nitrogen signal data.Area1

```

% All Nitrogen areas

arealCol = headers == "Areal";

nitroAreas = data(:,arealCol); % Nitroge signals (i.e. Areal)

nitroAreas = cell2num(nitroAreas);

% calibration (tin boats) blanks

calNitroBlankReps = nitroAreas(tinCalDex); % N blanks for calibration runs

calNitroBlank = nanmean(calNitroBlankReps); % Mean N blanks for calibration runs ONLY

allNitroCalBlanks(i) = calNitroBlank;

% sample tin circle blanks

sampleNitroBlankReps = nitroAreas(tinSampDex); % N tin blanks for sample runs

sampleNitroBlank = mean(sampleNitroBlankReps); % Mean N tin blanks for sample runs

allNitroSampBlanks(i) = nanmean(sampleNitroBlankReps);

% Acidified (Method) filter blanks

% filtNitroBlankReps = nitroAreas(acidFiltDex); % filter blanks
%
% filtNitroBlank = mean(filtNitroBlankReps);

```

Run CARBON calibration curves

```

% Nitrogen -- Amount of C in mg for sulfanilamide runs

weightCol = headers == "WeightVol";

weights = data(:,weightCol);

weights = cell2num(weights);

%      sulfCarbWeight  = (41.8/100).* weights(sulfDex); % mg C per sulfanilamide run

sulfNitroWeight  = (16.3/100).* weights(sulfDex); % mg N per sulfanilamide run

if ~isnan(calNitroBlank)

    sulfNitroAreas = nitroAreas(sulfDex) - calNitroBlank; % Sulfanilamide areas, blank corrected

else

    % Use mean cal blanks from previous runs if no blanks
    % available

    sulfNitroAreas = nitroAreas(sulfDex) - nanmean(allNitroCalBlanks(allNitroCalBlanks > 0)); % Sulfanilamide areas, blank corrected

end

if isnan(calNitroBlank) || isnan(sampleNitroBlank)
    warning('calNitroBlank and/or calNitroBlank were NaN; mean of previous runs used')
end

x = sulfNitroAreas;
y = sulfNitroWeight;

```

Warning: calNitroBlank and/or calNitroBlank were NaN; mean of previous runs used

Warning: calNitroBlank and/or calNitroBlank were NaN; mean of previous runs used

Calibration run QA

```

if strcmpi(calQA,'yes')

    areaNitroRatio = x./y; % area to nitroon ratio in sulfanilamide stds
    medianNitroRatio = nanmedian(areaNitroRatio);
    stdNitroRatio = std(areaNitroRatio);

```

```

absNitroResid = abs(areaNitroRatio - medianNitroRatio);
maxNitroResid = max(absNitroResid);
MADNitroRatio = mad(areaNitroRatio,1);
modZScores = 0.6745 * ((areaNitroRatio - medianNitroRatio) ./ MADNitroRatio);

% cal QA plot
% Tweaks so that figures renders the same in Mac & Linux

screens = handle(0); %
mainScreen = screens.MonitorPositions(1,3:4);

if ismac
    h0i = figure('Position',[mainScreen(1) * 1.06 378 364 636]);
elseif isunix
    h0i = figure('Position',[mainScreen(1) * 1.06 209 500 900]);
    h0i.Renderer = 'OpenGL';
end

subplot(2,1,1)

L1 = plot(areaNitroRatio,'o');
hold on

% axes properties

h1 = handle(gca);
h1.XLim = [0 ceil(numel(areaNitroRatio) * 1.06)];
h1.LineWidth = LineWidths;
grid on

h1.TickDir = 'out';
h1.TickLength = h1.TickLength .* 1.5;

L2 = hline(medianNitroRatio);
L2.LineStyle = '-';
L2.Color = newColors(2,:);

% An outlier is a value that is more than three SCALED median absolute ...
% deviations (MAD) away from the median:
%
% For a random variable vector A made up of N scalar observations,
% the median absolute deviation (MAD) is defined as

%     MAD' = median(abs(A-median(A)));
%
% The scaled MAD is defined as c * MAD
% where c = 1.4826 and is given by:
% -1/(sqrt(2)*erfcinv(3/2)).
%
% [***NOTE***] In Matlab mad(x) gives the MEAN absolute deviation from
% the media & while in R the function mad(x) gives the MEDIAN absolute
% deviation from the median. To get the scaled MAD in Matlab use the
% syntax mad(x,1)

% Graphical implementation of the above with the threshold values

X = [h1.XLim fliplr(h1.XLim)];
c = -1/(sqrt(2)*erfcinv(3/2)); %
M_i = scaledMADs * c * MADNitroRatio;

    y1 = medianNitroRatio - M_i;
    y2 = fliplr(medianNitroRatio + M_i);
Y = [y1 y1 y2 y2];
L3 = patch(X,Y,L2.Color);
L3.FaceAlpha = 0.3;
L3.EdgeColor = 'none';

% Y Lims

f = 4;
h1.YLim = [(medianNitroRatio - (M_i * f)) (medianNitroRatio + (M_i * f))];

%     L1 = plot(areaNitroRatio,'o');

%hline([meanNitroRatio meanNitroRatio-stdNitroRatio meanNitroRatio+stdNitroRatio],{'b', 'r', 'r'},{'mean+/-s.d.', '', ''})
xLab1 = xlabel('Run Order');
xLab1.Interpreter = 'latex';
yLab1 = ylabel('Signal:N ratio, $\frac{S_{\rm N}}{N}$ mg\, N');
yLab1.Interpreter = 'latex';
%title(['C calibration Q/A ' rundID ' run'])

% ID outlier Sulfanilamide runs (Mean )areaNitroRatio

goodCalRuns = abs(modZScores) < scaledMADs; % ~isoutlier(areaNitroRatio,'median',3); % index of cal runs to keep
badCalRuns = abs(modZScores) > scaledMADs; % isoutlier(areaNitroRatio,'median',3); % index of cal runs to keep

xx = 1:numel(areaNitroRatio);

L11 = plot(xx(badCalRuns),areaNitroRatio(badCalRuns),'o');
if ~isempty(L11)
    L11.Color = newColors(2,:);
end

x = x(goodCalRuns);
y = y(goodCalRuns);

```

```
% Legends

LegL3 = patch(h1.XLim(2) * [0.03 0.03 0.1 0.1],h1.YLim(1) + (range(h1.YLim) * [0.76 0.87 0.87 0.76]),L2.Color);
LegL3.FaceAlpha = 0.3;
LegL3.EdgeColor = 'none';
LegL4 = plot([h1.XLim(2) h1.XLim(2)] .* [0.03 0.1],h1.YLim(1) + (repmat(range(h1.YLim),1,2) .* [0.815 0.815]));
LegL4.LineWidth = LineWidths;
LegL4.Color = L2.Color;
ttl1 = text(h1.XLim(2) * 0.1167,h1.YLim(1) + (range(h1.YLim) * 0.815),'$x_{c-} < \hat{x} < x_{c+}$');
ttl1.Interpreter = 'Latex';
ttl1.FontSize = 18;


tt01 = text(h1.XLim(2) * 0.8,h1.YLim(1) + range(h1.YLim * 0.8),'a');
tt01.Units = 'Inches';
tt01.Position = tt01.Position;
tt01.FontSize = 20;

% Box
bx = patch(h1.XLim(2) * [0.018 0.018 0.55 0.55],h1.YLim(1) + (range(h1.YLim) * [0.73 0.9 0.9 0.73]),'w');
bx.LineWidth = LineWidths;


ti01 = title(runFileStr{1});
ti01.FontWeight = 'normal';
ti01.Interpreter = 'none';
ti01.FontSize = 12;

end
```


Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.



Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.



Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.



Plot Nitrogen Calibration curve

```
% C regression statistics

X = [ones(size(x,1),1) x];
[b,bint,r,rint,stats] = regress(y,X);
```

Linear regression of sulfNitro as a f(x) of C signal -- post QA

```
xx = union(0,x); % linear fit evaluation range

nitroModel = flipud(b)'; % polyfit(x,y,1); % Nitrogen mass vs C signal linear fit
%carbModel(2) = 0; % subtraction of intercept
z = polyval(nitroModel,xx);
```

```
subplot(2,1,2);

L33 = plot(x,y,'o');
L33.Color = newColors(1,:);

grid on
hold on

L4 = plot(sulfNitroAreas(goodCalRuns == 0),sulfNitroWeight(goodCalRuns == 0),'o');

if ~isempty(L4)
    L4.Color = newColors(2,:);
end

if sum(badCalRuns) > 0
    l0 = legend(L4,'Excluded','AutoUpdate','off');
    l0.Location = 'northwest';
end

aa = get(gca,'xlim');
bb = get(gca,'ylim');

L34 = plot(xx,z,'-');
L34.Color = newColors(2,:);
```

```

yL2 = ylabel('Nitrogen, mg');
yL2.Interpreter = 'Latex';

xL2 = xlabel('Nitrogen Signal, $S_{\rm C}$');
xL2.Interpreter = 'Latex';

% Display R^2 & regression eq on cal curve plot

tt01 = text(0.4 * aa(2), 0.22 * bb(2), ...
    ['$R^2 = ' num2str(stats(1),4) '$'], ...
    'fontSize', 16, ...
    'Interpreter', 'latex');
tt01.Units = 'Inches';
tt01.Position(1) = 1;
tt01.HorizontalAlignment = 'left';

newStr = regexp(num2str(b(2)), 'e-', 'split');

m = str2double(newStr{1});
m = round(m, 3);

newStr{1} = num2str(m);

regString = ['$ y = ' newStr{1} '\times 10^{-' newStr{2}(2) '}x + ' num2str(b(1),2) '$'];

tt02 = text(1.7, 0.17, ...
    regString, ...
    'fontSize', 16, ...
    'Interpreter', 'latex');
tt02.Units = 'Inches';
tt02.Position(1) = 1.0;
tt02.HorizontalAlignment = 'left';

h2 = handle(gca);
h2.LineWidth = LineWidths;
h2.TickDir = 'out';
h2.TickLength = h2.TickLength .* 1.5;


ttl = text(h2.XLim(2) * 0.8, h2.YLim(1) + range(h2.YLim * 0.75), 'b');
ttl.Units = 'Inches';
%ttl.Position = [3.3571 2.5904 0];
ttl.FontSize = 20;

% Plot position adjustment


%h2.Units = 'Inches';
h2.Position = h2.Position .* [1.25 1.4 0.95 1];

```


Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.




Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.



Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.



Warning: The DrawMode property will be removed in a future release. Use the SortMethod property instead.



Carbon and Nitrogen in filter samples

```

% Nitrogen

if ~isnan(sampleNitroBlank)
    sampleNitroAreas = nitroAreas(sampleDex) - sampleNitroBlank; % Sample Carbon areas blank corrected (minus tin square foil)
else
    sampleNitroAreas = nitroAreas(sampleDex) - mean(allNitroSampBlanks(allNitroSampBlanks > 0)); % Sample Carbon areas blank corrected (minus tin square foil)
end

%sampleCarbWeight = polyval(carbModel,sampleCarbAreas); % calculate sample Carbon in mg

% subtracting the intercept of cal curve
sampleNitroWeight = polyval(nitroModel,sampleNitroAreas) - nitroModel(2); % calculate sample Nitrogen in mg

% Nitrogen

% tile figures
if k == 1
    screenSize = get(0,'ScreenSize');

    plotWidth = 350;
    plotHeight = 325;

    nCols = floor(screenSize(3)/plotWidth);
    nRows = floor(screenSize(4)/plotHeight);

end

xPos = mod(k - 1,nCols) * plotWidth;
if k > nCols
    yPos = (screenSize(4) - (67 + ceil(k/nCols) * plotHeight));
else
    yPos = (screenSize(4) - (ceil(k/nCols) * plotHeight));
end

%

%% Export N concentration data to a file along with sample ID

```

```

        sampleID = str2double(runNames(sampleDex));

% if sampleID are not being converted to numbers try pass them directly

if sum(isnan(sampleID)) == numel(sampleID)
    sampleID = cell2num(runNames(sampleDex));

elseif any(isnan(sampleID))

    a = runNames(sampleDex);

    sampleID(isnan(sampleID)) = cell2num(a(isnan(sampleID)));
end

if k == 22
    stop
end

runNumber = repmat(k,numel(sampleID),1);

M = [runNumber sampleID sampleNitroWeight];

if k == 1
    dlmwrite('P06nitroWeighths.txt',M,...
        'delimiter',' ','roffset',0,'precision',9)

else

    dlmwrite('P06nitroWeighths.txt',M,'-append',...
        'delimiter',' ','roffset',0,'precision',9)
end

```

```

%% Export calibration curve and other performance metrics to file

runNumberStr = sprintf('runNumber_%d',k);

runDataPN.(runNumberStr).nitroModel = nitroModel;
runDataPN.(runNumberStr).calNitroBlankReps = calNitroBlankReps;
runDataPN.(runNumberStr).sampleNitroBlankReps = sampleNitroBlankReps;

if k == length(runFiles)
    save('P06PNrunDataStrct.mat','runDataPN')
end

```

```
end
```

Export figures to file

```

ag = findobj; % all graphical objects

nf = ag(2).Number; %max(ag==fix(ag)); % maximal integer value is number of figures

if strcmpi(exportFigs,'yes')

    cd([ homeStr 'Documents/CVO/Workshops/Writing/POCLatexVersion/Figures'])

    for i = 1:nf

        %print(i,'-dpdf',['Fig' num2str(i) 'PNCubeRun' rundID '.pdf'])
        print(i,'-dpng',['Fig' num2str(i) 'PNCubeRun' rundID '.png'])

    end

end

```

```

disp('Process completed')
toc

```

Process completed
Elapsed time is 226.891010 seconds.

Published with MATLAB® R2019b